



ArcelorMittal Bremen GmbH  
 Postfach 210220  
 28222 BREMEN  
 Carl-Benz-Strasse 30  
 28237 BREMEN  
 Telefon 0421/6480  
 Telefax 0421/6482251

A02

CERTIFICATE

Inspection certificate 3.1 EN 10204

A03

Page: 01 / 03

**20230029032-00**

A09DISPATCH NOTE

DATE SHIP  
 19.04.2023  
 DANICA HAV

A06.1 CUSTOMER

**FEON OY**  
 Teollisuuskatu 33  
 00510 Helsinki  
 Finland

A06.2 CONSIGNEE

**Feon OY**  
 Hepokarintie 9  
 23500 UUSIKAUPUNKI  
 Finland

0000325366

19.04.2023  
 DANICA HAV

A08.1 MANUFACTURER'S ORDER NR  
 2335398

DATE  
 27.01.2023

AGENCY'S ORDER NR.  
 GH32FEH003019

A07 CUSTOMER'S ORDER NUMBER

46004

A10 PART NUMBER

2463300

PRODUCT: hot rolled pickled coil

TERMS OF DELIVERY: EN 10051 (2010)


QUALITY: AMSTRONG®355MC

S355MC

STANDARD: AM FCE

EN 10149-2 (00-09-2013)

JITZ 1.82







ArcelorMittal Bremen GmbH  
Postfach 210220  
28222 BREMEN  
Carl-Benz-Strasse 30  
28237 BREMEN  
Telefon 0421/6480  
Telefax 0421/6482251

A02 CERTIFICATE  
Page: 03 / 03  
**20230029032-00**

Inspection certificate 3.1 EN 10204  
A05 ORIGINATOR OF THE DOCUMENT  
ArcelorMittal Bremen GmbH  
Carl-Benz-Straße 30, 28237 Bremen  
zeugnisse.bremen@arcelormittal.com  
Telefon 0421/6482813

A03

A09 DISPATCH NOTE  
DATE  
SHIP  
A08.1 MANUFACTURER'S ORDER NR  
DATE  
AGENCY'S ORDER NR.  
A07 CUSTOMER'S ORDER NUMBER  
A10 PART NUMBER

0000325366  
19.04.2023  
DANICA HAV  
27.01.2023  
GH32FEH003019  
46004  
2463300

A06.1 CUSTOMER  
FEON OY  
Teollisuuskatu 33  
00510 Helsinki  
Finland  
A06.2 CONSIGNEE  
Feon OY  
Hepokarintie 9  
23500 UUSIKAUPUNKI  
Finland

PRODUCT: hot rolled pickled coil  
TERMS OF DELIVERY: EN 10051 (2010)

QUALITY: AMSTRONG®355MC  
S355MC

STANDARD: AM FCE  
EN 10149-2 (00-09-2013)

JITZ 1.82

A08.2 ITEMTHICKNES	B09 mm	B10 WIDTH	B11 mm
19	5.00	1500.00	



QUALITY DEPARTMENT  
SITE EXPERT FOR INSPECTION  
Kochanowski

*Handwritten signature*

BREMEN 19.04.2023

We certify hereby that the delivery complies with the above mentioned specification.

C02 test direction relating to rolling direction (0°= L; 90°= T)	C04 specimen condition V:aged F:fresh N:normalised	C50 bend test 1:good	C43 result in J related to specimensize (C44)
--	--	----------------------	---